

FRET: Data analysis and protocol preparation

Chenyu Jin

Experimentalphysik I, Universitt Bayreuth, Bayreuth, Germany.
September 16, 2021

In the practice course, we are to measure the FRET efficiency in different cells with different methods. This is a guide for the data analysis and protocol preparation.

1 Sensitize Emission

In this part, we take for each cell three images:

- **Donor**, excited with donor excitation wavelength, and collect signal in donor emission channel;
- **s.e.** (sensitized emission), excited with donor excitation wavelength, and collect signal in acceptor emission channel;
- **Acceptor**, excited with acceptor excitation wavelength, and collect signal in acceptor emission channel;

Then we calculate FRET efficiency based on these three images and some correction factors.

This method, including the data analysis is not as straight-forward as it seems. There are many things to be pay attention to. Following is the explanation of some pitfalls and a recommended procedure of data analysis.

1.1 Dark current, background fluorescence and shading

Even if the signal-to-noise ratio on an image is high enough, we still need to clean up the signal before further analysis. In our case, what you will find out is, even if a ROI without any cell is




Figure 1: An example image of the cell, and the intensity profile along the yellow line, showing the non-zero background and the shading (the inclination of the base line).




Figure 2: Histogram of intensity when different ROIs are selected. This is an example from other construct that lights-up the whole cell. Please use only the bright membrane as the ROI in your analysis.

selected, the intensity measured by our detector is still not zero. This background is from (a) the dark current of the detector and (b) the background fluorescence on the slide. It adds to the real signal, and has to be substrated from the measured signal.

Another error is induced by the shading. Here it is mostly caused by the inhomogeneous illumination. In our case the effect brought by shading is minor so we will not correct shading in our analysis.

1.2 Selection of the region of interest (ROI)

Fluorescence is not homogeneous over the whole image. In our experiment, the membrane is the brightest, the cell body might also be fluorescent, and the surrounding should be dark. In a homogeneous region, the intensity distribution of pixels should be Gaussian (Poisson if too dark). Ideally, we should select only the “bright” homogeneous region as the ROI. Figure 2 shows what happens when some darker part is also included in the ROI.

Ideally, your ROI should be as homogeneous as possible, and also large. A Large ROI can average out the fluctuations, as you have already seen in the “Acceptor Photobleaching” experiment. **Therefore, please draw your ROI around the bright membrane part of the cells, and include a large area.**

1.3 Image processing protocol

The images can of course be processed with Python or Matlab. However, since we only need to process tens of images. I recommend Fiji (advanced version of ImageJ).

During the experiment, we have taken four images for each cell: **Donor**, **s.e.**, **Acceptor** and **BF** (bright field). They are not yet the **D**, **S** and **A** that we can feed to the equation to calculate the FRET efficiency. Some pre-processing is necessary:

- If to do image calculation, split frames (“Image – Stacks – stack to images”), and convert **Donor**, **s.e.**, **Acceptor** to 32-bits (“Image – Type”).
- Select a cell-free region on **BF**, measure the mean intensity on **Donor**, **s.e.**, **Acceptor**. These are the corresponding background values. You can use ROI manager to organize and save all the ROIs (“Image – Overlay”).
- If we only need a mean value, here we can draw ROI over the bright part, measure the mean intensity (“Analyze – Measure” or “M”, “Set Measurement” if something specific is needed), subtract the background values and feed the corrected intensity to the equations.
- To do image calculation, use “Process – Math” to substrate each image with the corresponding background value.




Figure 3: Steps of data processing, and the resulting histograms inside the red ROI. Before the image division, the **D** and **S** images are thresholded between [3, 250]. Please use only the bright membrane as the ROI in your analysis.

- The background subtraction will lead to some very small, even negative intensity values here and there on the image. You can hit “H” and see from the histogram. Later in image division these small values will cause problems, especially as denominator. We hence need a threshold to get rid of these small values. Open “Image – Adjust – Threshold”, set the lower threshold (upper bar) to around 2 or 3, click “Apply” and choose “Set to NaN”. Now all the pixels with a intensity value lower than the threshold will have a value “NaN”, which does not participate in calculation.

The result images are **D**, **S** and **A**. Use “Process – Math” or “Process – Image calculator” to calculate the “SE” and “FRET efficiency” images according to the equations. Choose a lookup table for the result FRET efficiency image (“Image – Lookup Tables”, e.g., “Red Hot”).

The FRET efficiency image now shows the FRET efficiency of each pixel. Please note that at each pixel it should be a number between 0 and 1. The whole image is actually a 2D plot. It is only meaningful when you add the calibration bar from “Analyze – Tools”. Don’t forget to add the scale bar as well.

Then please select on the FRET efficiency image the same ROI that you use to measure the mean intensity, and compare the mean FRET efficiency with what you calculated with the mean intensity.

1.4 What to include

In the first part, please include the following in your protocol:

- Tables with all the intensity and results of α , β , γ , and δ . Use the mean for further calculation.
- When calculating γ , please select the ROI once with only the bright membrane part, and once include the whole cell (with the dark inside and even some dark region outside). Also do calculation once without substrating the background. Then compare the standard deviation of all the results in this three cases.
- A table with all the results for SE and FRET efficiency. You can also estimate the mean distance of CFP and YFP in our cells.
- Pick 3 CY cells, do the image calculation to get the SE and FRET-efficiency images. Apply the same ROI used before on the image, and compare the mean of FRET-efficiency with the FRET-efficiency calculated from the mean intensity.

2 Donor emission after the acceptor bleaching

There is not much of data analysis to do in this part. In your protocol, please include:

- A typical CY image showing the ROIs.
- The corresponding plot of fluorescence intensity to time in donor and s.e. channel. Please pick similar colors for each ROI as shown on the image.
- Calculate the FRET efficiency with

$$E = 1 - \frac{D_{CY,pre}}{D_{CY,post}},$$

and list them in a table.

- Some description and discussion on the behaviour of the CFP and YFP cells.

3 Fluorescence Life-Time Measurement (FLIM)

In this part, we have measured the distribution of time between the excitation pulse and the emission photon. By fitting the histogram with exponential distribution we will get the life-time of the fluorescence.

There is only one excitation wavelength: 470 nm. Signals are collected in both donor (channel 2) and FRET (channel 1) channels. Hence for the three types of cells:

- C cells: donor signal in channel 2, leaked-through donor signal in channel 1.
- Y cells: cross-excited leak-through signal in channel 2, cross-excited signal in channel 1.
- CY cells: donor signal in channel 2, FRET signal in channel 1, mixed signal in channel 1+2. Note that only the dominating signals (should-be) are listed here.

The first task of data analysis is single exponential fit to the histograms (declining branch). The fitting result will be the life-time of the corresponding fluorescence. List all the results in tables.

The second task is double-exponential fit to the histogram of mixed signals (CY cells, channel 1+2): $f(t) = A_1 \cdot \exp(-t/\tau_1) + A_2 \cdot \exp(-t/\tau_2)$. There are 4 parameters to be found out with this fitting. Fitting with so many parameters will not give a meaningful result. The most famous example here is a quote by John von Neumann: “With four parameters I can fit an elephant, and with five I can make him wiggle his trunk”.

What we can do is to go back to the original purpose of fitting: to find out a reasonable function that is most close to the data. Here is one possible approach:

- Inherit A_1 from the fitting result of channel 1.
- Make A_2/A_1 the parameter to be found out. Take an initial guess from the fitting result of channel 2.
- Try a few numbers around the initial guess. Each time fit the data with only τ_1 and τ_2 with the fitting parameters. Calculate the abbreviation $\sum (y_{\text{data}} - y_{\text{fit}})^2$
- Find the A_2/A_1 at which the abbreviation is the smallest.

Do the double exponential fitting for one typical data set.

The third task is to fit the full data range. Also you only need to do it with one dataset, preferentially one from CY/channel 1, CY/channel 1+2, or Y/channel 1, because the signal-to-noise ratio is high.

The detected signal is the convolution (Deutsch: Faltung) of the real signal and the Instrument Response Function (IRF). So what you need is the convolution function (numpy.convolve, scipy.signal.convolve, or conv in Matlab). Due to the deadtime of the detector t_{tot} after the laser pulse, the measured signal $f(t)$ and the real signal $e(t)$ should be:

$$f(t) = \text{IRF} * e(t),$$

$$e(t) = \begin{cases} 0 & (t < t_{\text{tot}}) \\ A \cdot \exp(-t/\tau) & (t > t_{\text{tot}}) \end{cases}$$

A possible approach:

- Guess a t_{tot} from the measured signal (beginning of the rise, peak position, etc.).
- Take the result of the single exponential fit as the initial guess. You can fix A and aim to optimize τ .
- Try a few numbers around the initial guess. Calculate the abbreviation (as in the second task) after convolution, and find out the best fit.

Please do this with two types of IRFs: the one given by the TCSPC program, and a Gaussian curve deduced from the rising branch of the histogram. Please then plot the fitted $f(t)$, $e(t)$ together with the data set.

You are very encouraged to use other approaches, e.g., error function.

The FRET efficiency is then estimated with

$$E = 1 - \frac{\tau_{\text{donor in C}}}{\tau_{\text{donor in CY}}}$$

I know it doesn't work. Please add some discussion on why it fails.